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Instructions for Use

RIELLENS PROSTHETIC COMPONENTS FOR DENTAL IMPLANT

English

IFUA Rev01 of 10/2019

Analog

Riellens Analogs are of specific use to qualified professionals, with specific training. They are indicated for laboratory phase, used as replica of the implant platform.

1 - Intended Use

Riellens Analogs serve to reproduce the structure of the implant in order to assist prosthesis production or prosthetic planning in the plaster mold.

2 - Indication for Use

RIELLENS TABLE OF COMPONENTS			
COMPONENTS	INDICATION FOR USE	CONNECTION	DIAMETER
<p>Analog</p>	Used as replica, being a true copy of the platform and internal threads of the implant, it is attached to the plaster model for construction of the prosthesis that will be used on the patient.	HE	NP, RP, WP
		HI	Ø 3,3, 3,75, 3,4, 3,8, 4,1, 5,0
		TI	Ø 3,5/4,0, 4,5/5,0, NP, RP, WP
		CM	Ø2,9, 3,0, 3,3, 3,8, 4,5, NN, NP, RN, RP, WN, NC, RC

3 - Specifications and technical features

Riellens Analog are manufactured in Titanium or Stainless Steel in accordance with the standards ASTM F136, ISO 5832-3 (Standard Specification for Wrought Titanium-6Aluminum-4Vanadium) and ASTM F899 - 09 303 (Standard Specification for Wrought Stainless Steels for Surgical Instruments) ISO 5832-1 (Standard Specification for Implants for Surgery – Metallic materials Wrought Stainless Steels), have cylindrical format, are for single use and provided non-sterile.

4 - Form of Presentation

Packing

Riellens Analog is packed individually in blister packs with surgical grade paper, ensuring the integrity and cleaning of the product until the use by the professional. Information such as: description and size are available in the product's label.

CAUTION - Always use the package grip point, avoiding falls and direct contact with contaminated sites.

5 - Handling

After obtaining the mold, the analog must be screwed to the respective Transfer molding. To tighten the screw, use manual torque with digital key. After putting the plaster on the mold, covering the entire analog, move to the plaster placing phase, remove (separate) the plaster from the mold, thus obtaining the model. The analog reproduces exactly the position of the implant or prosthetic intermediary molded, serving as a basis for the prosthetic crown.

6 - Contraindications

As it is a material used in prosthetic laboratory, Riellens analog does not present any risk of complication to the patient.

7 - Components cleaning:

Riellens Analog are provided clean. The Analog used component that must be returned to the after-sales service, must be sent after cleaning without molding waste.

8 - Precautions and Warnings

- Note if the diameter is compatible with the prosthetic component.
- After casting, always evaluate the adaptation on the analog before cementing in the crown.
- Keep the integrity of the platform and the internal thread of the analog.
- Avoid excessive torque on the prosthetic screws.

WARNING
Dental product intended for use by licensed professionals. NON-STERILE.

9 - Storage Conditions

Always store at room temperature between 15° C and 35° C in its original packaging. Do not expose to direct sunlight.

10 - Product disposal

The legal procedures in force must be complied with for disposal of used or inadequate Analog and they must also have their characteristics removed by being cut, twisted or filed.

11 - Expiration Date and batch

Manufacture date, expiration and batch see package.

12 - Warning Instructions for Use

These INSTRUCTIONS FOR USE are available online at the manufacturer's website www.riellens.com.br/instrucao-de-uso/.

The INSTRUCTIONS FOR USE are indexed on the website by the TRADE NAME of the product, which can be found on the label of the product purchased. The INSTRUCTIONS FOR USE can be provided in printed material without additional cost if the professional prefers. Such request can be made at the Customer Service of the manufacturer.

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