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Instructions of Use RIELLENS PROSTHETIC COMPONENTS FOR DENTAL IMPLANT

English

IFUP Rev00 de 02/2017

Abutment/Mini Abutment

Riellens Abutments and Mini Abutments are for the specific use of qualified dentists, with specific training in implantology, including diagnosis, preoperative planning and surgical protocol.

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Riellens Abutments and Mini Abutments fix the prosthesis to the implant in order to reset or replace the dental structure and restore the patient's masticatory and aesthetic functions.

2. Indication for Use

Table 1 – Riellens Abutment/Mini Abutment Components

COMPONENTS	INDICATION FOR USE	MODELS	CONNECTION	DIAMETER	DIAMETER	HEIGHTS mm							
						1	2	3	4	5			
Mini Pilar Cônico	Used for making single or multiple bolted prosthesis being intermediate between the implant and the prosthesis	Straight	Used for making single bolted prosthesis or with more than one element (Multiple), being attached directly to the implant	HE	0 3,3, 4,1, 5,0								
				HI	0 3,4, 4,1, 5,0								
		Angled 17° and 30°	Used for making single bolted prosthesis or with more than one element (Multiple), with the need to correct the implant inclination with angulation of 17° and 30°, being attached directly to the implant	TI	0 3,5, 4,3, 5,0	X	X	X	X	X			
				CM	0 3,0, 0 3,3, 4,0, NN, NP, RN, RP, WN								
Pilar Cônico	Used for making single or multiple bolted prosthesis as intermediate between the implant and the prosthesis.	Straight	Used for making single bolted prosthesis or with more than one element (Multiple), being attached directly to the implant	HE	0 3,3, 4,1, 5,0								
				HI	0 3,4, 4,1, 5,0								
		Angled 17° and 30°	Used for making single bolted prosthesis or with more than one element (Multiple), with the need to correct the implant inclination with angulation of 17° and 30°, being attached directly to the implant	TI	0 3,5, 4,3, 5,0	X	X	X	X	X			
				CM	0 3,0, 0 3,3, 4,0, NN, NP, RN, RP, WN								
Pilar	Abutment Interface	Used for making single or multiple cemented zirconium-based prosthesis, specifically	Anti Rotational	Used for single prosthesis	HE	0 3,3, 4,1, 5,0							
					HI	0 3,4, 4,1, 5,0							
			Rotational	Used for prosthesis with more than one element (Multiple)	TI	0 3,5, 4,3, 5,0							
					CM	0 3,0, 0 3,3, 4,0, NN, NP, RN, RP, WN							
	Temporary Abutment	Used for making provisional prosthesis cemented, must be prepared in the laboratory before use.	Anti Rotational	Used for single prosthesis	HE	0 3,3, 4,1, 5,0							
					TI	0 3,5, 4,3, 5,0							
			Rotational	Used for prosthesis with more than one element (Multiple)	TI	0 3,5, 4,3, 5,0							
					CM	0 3,0, 0 3,3, 4,0, NN, NP, RN, RP, WN							
Definitive Straight Abutment	Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant.	Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant.	HI	0 3,4,4,1, 5,0									
					Definitive Angled Abutment	Used for making single and Multiple permanent cemented prosthesis, with the need to correct the implant inclination with angulation of 17°, being attached directly to the implant	TI	0 3,5, 4,3, 5,0	X	X	X	X	
							Immediate Temporary Abutment	Used for making temporary cemented prosthesis with more than one element (Multiple), and may be used immediately without the need of preparation. Single piece with external hexagonal fitting, where the digital key is introduced to give the torque.	CM	0 3,0, 0 3,3, 4,0, NN, NP, RN, RP, WN			
					Definitive Solid Abutment	Single body part. Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant			CM	RN, WN			

3. Technical specifications and characteristics

Riellens Abutments and Mini Abutments are manufactured in titanium according to standards ASTM F136 and ISO 5832-3 (Standard Specification for Wrought Titanium-6Aluminum-4Vanadium), intended for single use, provided non-sterile and have a cylindrical in shape.

4. Method of Presentation

Package

Riellens Abutments and Mini Abutments are individually packaged in blister packaging, ensuring the integrity and cleanliness of the product until the use by the professional. Information such as description, dimension and occasionally angulation are available on the product label. The components Conical Abutment, Mini Conical Abutment, Abutment Interface, Temporary Abutment, Immediate Temporary Abutment, Definitive Straight Abutment, Definitive Angled Abutment, Definitive Solid Abutment, Definitive Universal Abutment, Definitive Rulator Abutment, Definitive Cementable Abutment and Definitive Bar and Bridge Abutment are accompanied by three labels of traceability, which must mandatorily be fixed: (i) in the medical record; (ii) in the document to be delivered to the patient; (iii) the invoice of product

ATTENTION – Always use the handle spot of the package, avoiding falls and direct contact with contaminated places.

In the event of damage to the package (blister) or after the expiration date, Riellens Abutments and Mini Abutments cannot be used. All Abutments and Mini Abutments are sold with their respective screw, which can be sold separately as a spare part.

5. Recommended Torques/Used Wrenches

Table 2 – Correspondence Table Abutment/Mini Abutment x Wrench

PROSTHETIC COMPONENTS	INDICATION FOR USE	TYPES OF THREADS	DIAMETER (Ø)	RECOMMENDED TORQUE	MAXIMUM ALLOWABLE TORQUE	USED WRENCH		
						HEX 1.2	UNIG (Tor)	SIZOR RE 1.3
	Indicated for fixing components directly to the implant, with controlled torque application providing a perfect adaptation. They are made with Titanium ASTM F136 in single body having at one end an external thread and at the other a Hexagonal, Square or Unigrip inner socket.	Thread M1.4	RP, 3.0	20 N.cm	26 N.cm	X	X	-
		Thread M1.6	3.3, 3.4, 4.0, 4.1, 5.0,	25 N.cm	32.5 N.cm	X	X	-
		Thread M1.8, M2.0 y M2.5	WP, 3.4, 3.5, 4.3, 5.0	32 N.cm	41.6 N.cm	X	-	X
	Indicated for fixing components directly to the implant, with controlled torque application providing a perfect adaptation. They are made with Titanium ASTM F136 in single body having at one end an external thread and at the other a Hexagonal, Square or Unigrip inner socket.	Specific for Mini Conical Abutment and Conical Abutment with Thread M1.6	3.3, 4.1	20 N.cm	26 N.cm	x	-	-
		Specific for Mini Conical Abutment with Thread M1.8 M2.0	3.5, 4.1	25 N.cm	32.5 N.cm	-	X	-
	Indicated for fixing components directly to the implant, with controlled torque application providing a perfect adaptation. They are made with Titanium ASTM F136 in single body having at one end an external thread and at the other a Hexagonal, Square or Unigrip inner socket.	Specific for Mini Conical Abutment with Thread M2.5	5.0	25 N.cm	32.5 N.cm	-	X	-
		Specific for Conical Abutment with Thread M1.8 M2.0	3.5, 4.1, 4.3	25 N.cm	32.5 N.cm	X	-	-
	Used for making temporary cemented prosthesis with more than one element (Multiple), and may be used immediately without the need of preparation. Single piece with external hexagonal fitting, where the digital key is introduced to give the torque.	Thread M1.6 y M1.8	3.3, 3.4, 3.5, 4.1 y 5.0 NP	25 N.cm	32.5 N.cm	X	-	-
		Thread M2.0 y M2.5	3.4, 4.1, 4.3 y 5.0 RP, WP, RN, WN	32 N.cm	41.6 N.cm	X	-	-

6. Handling

Place the Abutment or Mini Abutment over the implant and apply the recommended torque (specific to each prosthetic interface). To cement the metal Abutments with the prosthesis, use the cement according to the quantity indicated by the manufacturer:
Definitive Cementation: Zinc phosphate.
Temporary Cementation: Calcium hydroxide

7. Contraindications

The use of Riellens Abutments and Mini Abutments is contraindicated for patients who present chronic periodontal inflammation, inadequate oral hygiene, intra-oral infection, insufficient bone quality, peri-implant infection, endodontic or periodontal problems, untreated periodontitis, Tobaccoism and allergy to chemical elements comprising its raw material, such

8. Cleaning of the Components

The use of Rielens Abutments and Mini Abutments is contraindicated for patients who present chronic periodontal inflammation, inadequate oral hygiene, intra-oral infection, insufficient bone quality, peri-implant infection, endodontic or periodontal problems, untreated periodontitis, Tobaccoism and allergy to chemical elements that make its raw material, such as V, Ti, Al, O, N, Fe, H and C. Rielens provides all its Abutments and Mini Abutments sanitized. Any component used that must be returned to the after-sales service shall be shipped sterile after pre-disinfection, cleaning and decontamination in the office.

9. Sterilization

The components Conical Abutment, Mini Conical Abutment, Abutment Interface, Temporary Abutment, Immediate Temporary Abutment, Definitive Straight Abutment, Definitive Angled Abutment, Definitive Solid Abutment, Definitive Universal Abutment, Definitive Rilatator Abutment, Definitive Cementable Abutment and Definitive Bar and Bridge Abutment are non-sterile, thus they shall be sterilized in an autoclave, according to the professional protocol and manufacturer's instruction before being used on the patient. Rielens recommends placing the components one by one in a self-adhesive envelope for steam sterilization in an autoclave, moist heat at 134°C +/- 1°C for 18 minutes, resterilization is forbidden.

ATTENTION
Dental product intended for the use of qualified professionals. Supplied **NON-STERILE** and must be sterilized before use.

10. Precautions and Warnings

Failure to follow the procedures outlined in this instruction for use may lead to one or all of the following complications: implant or component damage, components loosening, inadequate final restoration or poor adaptation of mobile and fixed prostheses, implant failure. Rielens RECOMMENDS SINGLE USE. RESTERILIZATION IS FORBIDDEN. The reuse of the components will result in damage to the patient's health who may come in contact with a contaminated product.

The use of components or instruments that are not compatible will render all warranties void. The correct use and handling of this product is the sole responsibility of the user. Rielens takes no responsibility for damages resulting from incorrect use of its products. The implementation of the prosthesis must be done in a surgical environment, which must be clean and suitable for the manipulation and use of health products.

11. Warnings

Torsion superior to the recommended torque may result in damage to the screw fitting and fracture, causing prosthesis failure.
Torsion inferior to the recommended torque may result in component loosening, which may lead to prosthesis failure.

12. Attention

Use the calibrated torque meter until the recommended torque is obtained to tighten the screw and use the fixing wrench, according to table 2, to stabilize. Follow the Torquemeter manufacturer's instructions.
Ensure that the connection screws between the implant and the abutment are correctly aligned on all secondary components to prevent wear on the threads.

13. Adverse Effects

Possible Complications:
- Pain, swelling, phonetic difficulties, gingival inflammation, localized or systemic infections, aesthetic problems, restoration or dental fractures.

14. Post-implementation Care

Advise the patient to see the dental surgeon periodically.

15. Condition of Storage

Always store at room temperature between 15°C and 35°C in the original package. Do not expose to direct sunlight.

16. Disposal of the Product

For disposal of Abutments and Mini Abutments used or considered inadequate, it is necessary to follow the legal sanitary procedures in force for the disposal of solid wastes in the

17. MRI - Magnetic Resonance Safety

The non-clinical test demonstrated that the product is MR conditional. A patient with this device may be safely checked on an MR system that meets the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla only.

- Magnetic wires of maximum spatial gradient of 4000 Gauss/cm or minimum of (40T/m).

- Maximum reported MR system, full body average specific absorption rate (SAR) of 2W/kg (Normal Operating Mode) or 4W/kg (First Level Controlled Mode).

Under the scan conditions defined above, the product is expected to produce a maximum temperature rise of 4.1°C after 15 minutes of continuous operation. In non-clinical tests, the image artifact caused by the device extends approximately 30mm from the product when represented with a sequence of gradient echo pulses with a magnetic resonance system of 3 Tesla.

If there is no MR symbol on the product label, note that the product has not been evaluated for safety and compatibility in the MR environment. The product has not been tested for heating or migration in the MR environment.

18. Expiration Date and Batch

See the manufacture date, expiration and batch on the package.

19. INSTRUCTION FOR USE WARNING

This INSTRUCTION FOR USE is available in non-printed (electronic) format, through the electronic address of the manufacturer at www.riellens.com.br. The INSTRUCTIONS FOR USE are indexed on the website through the respective product TRADE NAME, informed on the label of the product purchased.

If it is in the professional's interest, the INSTRUCTIONS FOR USE may be provided in printed form, at no additional cost. We ask the professionals to notify and inform Rielens of any problems that may occur with the product during its use. The requests must be made to the manufacturer's SAC (Customer Service – Telephone 55 11 – 4158-5018).

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