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Instructions for Use

RIELLENS PROSTHETIC COMPONENTS FOR DENTAL IMPLANT

English

IFUT Rev01 10/2019

Transfer

Riellens Transfers are used specifically by qualified professionals, with specific training in implant dentistry, including diagnosis, pre-surgery planning and surgical protocol.

1. Intended Use

The Riellens Transfers serve as transfer (molding) on Conical Abutment, Mini Conical Abutment and Implants. It is indicated to open molder (drag) or closed molder (spare) techniques.

2. Indication for Use

COMPONENTS TRANSFERS RIELLEN'S					
COMPONENTS	INDICATION FOR USE	MODELS	FEATURES	CONNECTION	DIAMETER
 Transfer in Stainless Steel	Used in the transfer of the dental implant placement installed in the oral cavity to the plaster model, through a molding.	Open Tray	Used for single and multiple prosthetic moldings with open tray	HE	Ø NP, RP, WP
		Closed Tray	Used for single and multiple prosthetic moldings with Closed tray	HI	Ø 3.75
 Transfer in titanium				TI	Ø 3.3
				CM	Ø 3.0, Ø 3.8, Ø 4.5, Ø 4.5/5.0, NN
		Open Tray	Used for single and multiple prosthetic moldings with open tray	HI	Ø 3.8, 7.3, NP, RP, WP
		Closed Tray	Used for single and multiple prosthetic moldings with Closed tray	TI	Ø 3.5/4.4, 4.5/5.0, NP, RP, WP
				CM	Ø NN, NP, RN, RP, WN, NC, RP

3. Specifications and technical features

Riellens Transfers are manufactured in titanium in accordance with the standards ASTM F136 and ISO 5832-3 (Standard Specification for Wrought Titanium-6Aluminum-4Vanadium) or Stainless Steel ASTM F899-09 303 (Standard Specification for Wrought Stainless Steels for Surgical Instruments) ISO 5832-1 (Standard Specification for Implants for Surgery – Metallic materials Wrought Stainless Steels), have cylindrical format, are reusable and provided non-sterile.

They have two parts, Transfer + Screw and are available as open molder and closed molder in relation to the molding technique.

4. Form of Presentation

Packing

Riellens Transfers are packed individually in blister packs, ensuring the integrity and cleaning of the product until the use by the professional. Information such as: description, size, manufacturing date and batch are available in the product's label.

CAUTION - Always use the package grip point, avoiding falls and direct contact with contaminated sites.

Riellens Transfers may not be used if the package (Blister) is damaged or after the expiration date.

All Riellens Transfers are marketed with their respective screw, and may be sold separately as spare part.

5. Handling

Place the Transfer selected according to the planned molding technique and install directly on the implant or intermediary. To tighten the screw, use manual torque with digital key.

6. Contraindications:

This product does not present contraindications, as long as used correctly for the purposes indicated.

7. Components cleaning:

Riellens Transfers are provided clean. Any used component that must be returned to the after-sales service must be sent after cleaning and decontamination in the dental clinic.

The components must be cleaned, disinfected and sterilized before each use.

- Put the pieces in the cleaning bath for removal of organic matter from the components as soon as possible after use; Use enzymatic detergent, following the instructions for dilution and exposure time determined by the manufacturer.

- Perform the mechanical removal of the filth of the components by manual brushing (soft brushes) or vibrations produced by ultrasonic cleaner.

- Wash with clean water to remove chemical waste and detergents.

- Dry the component always using disposable soft lint-free cloth.

Note - The responsible must be equipped with individual protection equipment (gloves, masks, goggles, aprons, hats, etc.).

8. Sterilization

The Transfers and Screws components are not sterile, thus they should be autoclaved as professional protocol and manufacturer's instructions before reuse on the patient.

Riellens recommends placing the components in the proper package for sterilization in autoclave. According to the following table:

Method	Conditions
Autoclave, Humid heat	134° for 18min

9. Precautions and Warnings
 The use of incompatible components and instruments voids all warranties. Users are responsible for the proper use and handling of this product. Riellens will not be responsible for any damage caused due to inadequate use of its products.
 The prosthesis installation environment must be clean and appropriate to the handling and use of health products.

WARNING
Dental product intended for use by licensed professionals. Provided NON-STERILE.
Must be sterilized before each use.

10. Caution
 To tighten the screw use a digital manual key until the desired torque.
 Ensure the transfers screws are properly aligned between the implant/and or the intermediate to prevent wear on the threads.

11. Adverse Effects
 This product does not present adverse events, as long as used correctly for the purposes indicated.

12. Storage Conditions
 Always store it at room temperature between 15°C and 35°C in the original package. Do not expose to direct sunlight.

13. Product disposal
 The legal procedures in force must be complied with for disposal of used or inadequate Transfers and they must also have their characteristics removed by being cut, twisted or filed.

14. Expiration Date and batch
 Manufacture date, expiration and batch see package. This product shall be disposed if it loses its functionality.

15. WARNING INSTRUCTIONS FOR USE
 These INSTRUCTIONS FOR USE are available online at the manufacturer's website www.riellens.com.br
 The INSTRUCTIONS FOR USE are indexed on the website by the TRADE NAME of the product, which can be found on the label of the product purchased.
 The INSTRUCTIONS FOR USE can be provided in printed material without additional cost if the professional prefers. Such request can be made at the Customer Service of the manufacturer.

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