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Instruction for Use RIELLEN'S PROSTHETIC COMPONENTS FOR DENTAL IMPLANT

English

FUC Rev02 of 10/2020

Heling/Protective Cover

Riellens Helings and Protective Covers are for the specific use of qualified professionals, with specific training in implantology, including diagnosis, preoperative planning and surgical protocol.

1. Intentions of Use

Helings and Protective Covers have the function of exposing the implant to the oral environment, preserving and increasing, if necessary, the inserted gum and adjust the thickness of soft tissues, healing the peri-implant gingival tissue.

2. Table 1 - Indication for Use

RIELLEN'S HEALING COMPONENTS/PROTECTIVE COVER										
COMPONENTS	INDICATION FOR USE	CHARACTERISTICS	CONNECTI ON	DIAMETER	HEIGHTS mm					
					2	3	4	5	6	7
	The Heling and Protective Cover are used in the second phase of implant surgery, with the purpose of creating a cylindrical space in the gum, performing appropriate healing of the gingival tissue.	It is a single body component applied directly on the implant	HE	NP, RP, WP	x	x	x	x	x	x
			HI	NP, RP, WP, Ø 3,75, 4,1, 4,3, 5,0, 6,0	x	x	x	x	x	x
			TI	NP, RP, WP, 6,0	x	x	x	x	x	x
	It is a single body component applied directly on the intermediate part	It is a single body component applied directly on the intermediate part	HE	NP, RP, WP	x	x	x	x	x	
			HI	NP, RP, WP, Ø 3,75, 6,0	x	x	x	x	x	
			CM	Ø 2,9, 3,3, 3,8, 4,0, 4,5, RN, WN, NC, RC, NP,	x	x	x	x	x	x

KEY: HE: External Hexagon; HI: Internal Hexagon; CM: Morse Cone; TI: Inner Triangle



3. Technical specifications and characteristics

Helings and Protective Covers are made of titanium according to ASTM F136 and ISO 5832-3 (Standard Specification for Wrought Titanium-6Aluminum- 4Vanadium), they have a cylindrical shape and present a socket shaft (1.2 mm) or Unigrip at the upper end, and are provided non-sterile.

4. Way of Presentation

Packaging

Riellens Helings are individually packaged in blister packaging, ensuring the integrity and cleanliness of the product until the use by the professional. Information such as description and fitting are available on the product label.

5. Handling

The Heling and protective cover must be used for a short term, and it is recommended to remove them within 30 days after implantation. For greater safety, protect surgical drape, avoiding that if the product is released into the mouth, it will not be swallowed or aspirated by the patient.

Heling - Place the Heling over the implant and apply recommended torque (specific to each prosthetic interface).

Protective cover - Place the protective cover over the mini abutment or conical abutment and apply recommended torque (specific to each prosthetic interface).

6. Recommended Torques/Used Wrenches

Table 2 - Correspondence Table Heling/Protective Cover x Wrench

COMPONENTS	TYPES OF THREADS	PLATFORM	RECOMMENDED TORQUE	MAXIMUM ALLOWABLE TORQUE	WRENCH USED	
					HEX. 1.2	UNIG (Torx)
	M1.4 Thread	HE CM MU (RP) BL (2.9)	20 N.cm	26 N.cm	X	X
	M1.6 Thread	HE HI CM HI CE (NP, RP, WP, 6.0) BL (NC, RC); NA (NP) GMN (Ø 3.3, 4.5) BLR (3.8)	25 N.cm	32.5 N.cm	X	X
	M1.8 Thread	HE HI TI CM MU (WP) MS (Ø3,75) RS (NP) CMN (Ø3.3, 4.5) CMD (Ø3.8) MGA (4,0)	25 N.cm	32.5 N.cm	X	X
	M2.0 Thread	HE TI CM 3HE (NP, WP), HEB (RP) RS (RP, WP, 6.0) TL, ST, SI (RN, WN) NA (RP) CMD (Ø4.5/5.0)	32 N.cm	41.6 N.cm	X	X
	M2.5 Thread	HE TI HEB (WP) RS (RP, WP)	32 N.cm	41.6 N.cm	X	X

7. Contraindications:

The use of Riellens Helings and Protective Covers are contraindicated in cases of patients who are allergic to any chemical element that make the raw material such as: V, Ti, Al, O, N, Fe, H and C.

8. Adverse Effects/Technical Complaints

Possible complications: Gingival inflammation, localized infections, maladaptation, mechanical problems

9. Component cleaning:

Riellens makes available all its Helings and Protective covers sanitized.

10. Sterilization

These products are provided non-sterile; thus, they shall be sterilized in an autoclave.

- Use chemical and biological indicators for each sterilization cycle, according to practice protocol.
- Allow Riellens Helings/Protective covers to dry and cool in the autoclave before handling to avoid contamination.
- Riellens recommends sterilization at the parameters of 134 °C for 15 minutes.
- Riellens recommends placing the components in their own packaging for sterilization in an autoclave.

11. Precautions and Warnings

- The use of components or instruments that are not compatible will render all warranties void. The correct use and handling of this product is the sole responsibility of the user. Riellens takes no responsibility for damages resulting from incorrect use of its products.
- The environment for installing the prosthesis must be clean and suitable for the handling and use of health products.
- IT IS "FORBIDDEN TO REPROCESS" the Heling and Protective cover; the reuse of the components will result in damage to the patient's health who may come in contact with a contaminated product.
- Torque superior to the recommendation may result in damage to the Heling and Protective covers fitting and fracture, causing prosthesis failure. Torque inferior to the recommendation may result in component loosening, which may lead to poor healing of the gingival tissue.

ATTENTION
Dental product intended for use by qualified professionals.
Provided NON-STERILE, never reuse Riellens Heling/Protective cover risk of infection with infectious agents.
Caution when handling the component; risk of accidentally aspirating or swallowing.

12. Alerts / Cautions

- In case of technical complaints or adverse events with Riellens products, please notify us through Customer Service.
- Ensure that the Heling and Protective covers are correctly aligned to prevent thread wear.
- In the event of damage to the package (blister) or after the expiration date, Riellens Heling and Protective covers cannot be used.
- Any used components that must be returned to the after-sales service should be shipped sterile after pre-disinfection, cleaning and decontamination in the practice.

13. Storage Conditions

Always store at room temperature between 15 °C and 35 °C in the original package. Do not expose to direct sunlight.

14. Product disposal

For the disposal of the Heling and Protective covers used or considered inadequate, it is necessary to follow the current legal procedures for the disposal of products, and they must be defused by being cut, bent or filed away.

15. Expiration Date and Batch

See the manufacture date, expiration and batch on the package.

16. INSTRUCTION FOR USE ALERT

This INSTRUCTION FOR USE is available in non-printed (electronic) format, through the electronic address of the manufacturer at www.riellens.com.br. The INSTRUCTIONS FOR USE are indexed on the website through the respective product TRADE NAME, informed on the label of the product purchased. If it is in the professional's interest, the INSTRUCTIONS FOR USE may be provided in printed form, at no additional cost. The request must be made to the SAC (Customer Service) of the manufacturer.

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