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# Instructions of Use

## RIELLENS PROSTHETIC COMPONENTS FOR DENTAL IMPLANT

English

IFUP Rev01 de 09/2019

### Abutment/Mini Abutment

Riellens Abutments and Mini Abutments are for the specific use of qualified dentists, with specific training in implantology, including diagnosis, preoperative planning and surgical protocol.

#### 1. Intended Use

Riellens Abutments and Mini Abutments fix the prosthesis partial or total to the implant in order to reset or replace the dental structure and restore the patient's masticatory and aesthetic functions.

#### 2. Indication for Use

Table 1 – Riellens Abutment/Mini Abutment Components

COMPONENTS	INDICATION FOR USE	MODELS	CHARACTERISTICS	COMPATIBLE WITH IMPLANT SYSTEM	DIAMETER	HEIGHTS mm					
						1	2	3	4	5	
Mini Conical Abutment	Used for making single or multiple bolted prosthesis being intermediate between the implant and the prosthesis	Mini Conical Abutment	Used for making single bolted prosthesis or with more than one element (Multiple), being attached directly to the implant	HE	Ø NP, RP, WP						
				HI	Ø NP, RP, WP	X	X	X	X	X	
				TI	Ø 3.75, NP, RP, WP						
Mini Conical Angled Abutment 17° and 30°	Used for making single bolted prosthesis or with more than one element (Multiple), with the need to Correct the implant inclination with angulation of 17° and 30°, being attached directly to the implant	Mini Conical Angled Abutment 17° and 30°	Used for making single bolted prosthesis or with more than one element (Multiple), with the need to Correct the implant inclination with angulation of 17° and 30°, being attached directly to the implant	CM	Ø 3.5, 3.8, 4.5 NN, NP, RN, RP, WN						
				HE	Ø NP, RP, WP						
				HI	Ø NP, RP, WP						
Conical Abutment	Used for making single or multiple bolted prosthesis as intermediate between the implant and the prosthesis.	Conical Abutment	Used for making single bolted prosthesis or with more than one element (Multiple), being attached directly to the implant.	HE	Ø NP, RP, WP						
				HI	Ø NP, RP, WP	X	X	X	X	X	
				TI	Ø 3.75, NP, RP, WP						
Conical Angled Abutment 17° and 30°	Used for making single bolted prosthesis or with more than one element (Multiple), with the need to Correct the implant inclination with angulation of 17° and 30°, being attached directly to the implant.	Conical Angled Abutment 17° and 30°	Used for making single bolted prosthesis or with more than one element (Multiple), with the need to Correct the implant inclination with angulation of 17° and 30°, being attached directly to the implant.	CM	Ø 3.5, 3.8, 4.5 NN, NP, RN, RP, WN						
				HE	Ø NP, RP, WP						
				HI	Ø NP, RP, WP						
Abutment	Interface Abutment	Interface Abutment	Used in the manufacture of single or multiple cemented prosthesis specific for CAD / CAM system being fixed directly to the implant.	Anti Rotational	Used for single prosthesis	HE	Ø NP, RP, WP				
				Rotational	Used for prosthesis with more than one element (Multiple)	HI	Ø 3.3, 3.7, 3.75, 4.2, 4.5, NP, RP, WP.				
				CM	Ø 3.3, NP, RP, RP.						
Temporary Abutment	The abutment is connected directly to the dental implant endosseous, used in the manufacture of prosthesis for cemented implant unitary or multiple	Temporary Abutment	Used in the manufacture of temporary cemented prosthesis, it must be prepared in the laboratory before use.	Anti Rotational	Used for single prosthesis						
				Rotational	Used for prosthesis with more than one element (Multiple)						
				HE	Ø NP, RP, WP						
Definitive Straight Abutment	Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant.	Definitive Straight Abutment	Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant.	HE	Ø NP, RP, WP						
				HI	Ø 4.1, 4.3, 5.0, NP, RP, WP						
				TI	Ø 3.3, NP, RP						
Definitive Angled Abutment	Used for making single and Multiple permanent cemented prosthesis, with the need to correct the implant inclination with angulation of 17°, being attached directly to the implant.	Definitive Angled Abutment	Used for making single and Multiple permanent cemented prosthesis, with the need to correct the implant inclination with angulation of 17°, being attached directly to the implant.	CM	Ø 3.0, NN, NP, RN, RP, WN	X	X	X	X	X	
				HE	Ø NP, RP, WP						
				HI	Ø 4.1, 4.3, 5.0, NP, RP, WP						
Immediate Temporary Abutment	Used for making temporary cemented prosthesis with more than one element (Multiple), and may be used immediately without the need of preparation. Single piece with external hexagonal fitting, where the digital key is introduced to give the torque.	Immediate Temporary Abutment	Used for making temporary cemented prosthesis with more than one element (Multiple), and may be used immediately without the need of preparation. Single piece with external hexagonal fitting, where the digital key is introduced to give the torque.	HE	Ø NP, RP, WP						
				HI	Ø 4.1, 4.3, 5.0, NP, RP, WP						
				TI	Ø 3.3, NP, RP						
Definitive Solid Abutment	Single body part. Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant	Definitive Solid Abutment	Single body part. Used for making permanent prosthesis cemented with more than one element (Multiple), being attached directly to the implant	CM	RN, WN				X	X	
				HE	Ø NP, RP, WP						
				HI	Ø 4.1, 4.3, 5.0, NP, RP, WP						

LEGEND: HE: External Hexagon; HI: Internal Hexagon; CM: Cone Morse; IT: Internal Triangle

WARNING: All Riellens Pillars and Mini Pillars have dimensionally compatible fittings with the HE, HI, CM and TI connection implants on the market.



#### 3. Technical specifications and characteristics

Riellens Abutments and Mini Abutments are manufactured in titanium according to standards ASTM F136 and ISO 5832-3 (Standard Specification for Wrought Titanium-6Aluminum-4Vanadium), intended for single use, provided non-sterile and have a cylindrical in shape.

#### 4. Method of Presentation

##### Packing

Riellens Abutments and Mini Abutments are individually packaged in blister packs ensuring product integrity and cleanliness until professional use. Information such as description, size and possibly angulation are available on the product label. Components Conical Abutment, Mini Conical Abutment, Interface Abutment / Tibase Temporary Abutment, Immediate Temporary Abutment, Abutment; Solid Pillar, are accompanied by three traceability labels, which must be fixed: (i) in the medical record; (ii) the document to be delivered to the patient; (iii) in the purchase invoice of the product.

In the event of damage to the primary packaging (blister) or after the expiration date of the Abutments and Mini Abutments, Riellens does not guarantee product cleanliness

All Abutments and Mini Abutments are sold with their respective screw, which can be sold separately as a spare part.

ATTENTION – Always use the handle spot of the package, avoiding falls and direct contact with contaminated places.

#### 5. Recommended Torques/Used Wrenches

Table 2 – Correspondence Table Abutment/Mini Abutment x Wrench

PROSTHETIC COMPONENTS	INDICATION FOR USE	TYPES OF THREADS	DIAMETER (Ø)	RECOMMENDED TORQUE	MAXIMUM ALLOWABLE TORQUE	USED WRENCH			
						HEX 1.2	UNIG (Torx)	SQUARE 1.3	
	Indicated for fixing components directly to the implant, with controlled torque application providing a perfect adaptation. They are made with Titanium ASTM F136 in single body having at one end an external thread and at the other a Hexagonal, Square or Unigrip inner socket. . .	Thread M1.4	RP, 3.0	20 N.cm	26 N.cm	X	X	-	
		Thread M1.6	3.3, 3.4, 4.0, 4.1, 5.0,	25 N.cm	32.5 N.cm	X	X	-	
		Thread M1.8, M2.0 y M2.5	RP, NP, WP, 3.4, 3.5, 4.3, 5.0, RN, WN	32 N.cm	41.6 N.cm	X	-	X	
	Indicated for fixing components directly to the implant, with controlled torque application providing a perfect adaptation. They are made with Titanium ASTM F136 in single body having at one end an external thread and at the other a Hexagonal, Square or Unigrip inner socket. . .	Specific for Mini Conical Abutment and Conical Abutment with thread M1.6	3.3, 4.1, 5.0	20 N.cm	26 N.cm	x	-	-	
			Specific for Mini Conical Abutment with thread M1.8 M2.0	3.5, 4.1, 4.3, 5.0	25 N.cm	32.5 N.cm	X	-	-
				5.0			-	X	-
	Indicated for fixing components directly to the implant, with controlled torque application providing a perfect adaptation. They are made with Titanium ASTM F136 in single body having at one end an external thread and at the other a Hexagonal, Square or Unigrip inner socket. . .	Specific for Conical Abutment with thread M1.8 M2.0	3.3, 3.4, 3.5, 4.1, 4.3, 5.0	25 N.cm	32.5 N.cm	X	-	-	
			3.4, 3.5, 4.1, 4.3, 5.0			X	-	-	
	Used for making temporary cemented prosthesis with more than one element (Multiple), and may be used immediately without the need of preparation. Single piece with external hexagonal fitting, where the digital key is introduced to give the torque.	Thread M1.6 y M1.8	3.3, 3.4, 3.5, 4.1 and 5.0 NP	25 N.cm	32.5 N.cm	X	-	-	
			Thread M2.0 y M2.5			3.4, 4.1, 4.3 and 5.0 RP, WP, RN, WN	32 N.cm	41.6 N.cm	X

#### 6. Handling

Clinical procedure after prosthetic process:

Screw Prosthesis:

1. Clean and disinfect the final restoration as directed by the restorative material manufacturer.
2. Remove the healing or temporary restoration from the patient.
3. Attach and tighten the permanent reset with a digital key. It is recommended to check the adaptation of the abutment using radiographic imaging.
4. Tighten the Abutment with the torque wrench, applying the recommended torque as shown in Table 2 - Matching Abutment / Mini Abutment x Wrench

Cemented Prosthesis:

1. Clean and disinfect the final restoration as directed by the restorative material manufacturer.
2. Remove the healing or temporary restoration from the patient.
3. To cement the metal abutments with the prosthesis, use the cement according to the amount indicated by the manufacturer:
  - Ultimate Cmentation: Zinc phosphate.
  - Provisional Cmentation: Calcium hydroxide

## 7. Contraindications

The use of Rielens Abutments and Mini Abutments is contraindicated for patients who present chronic periodontal inflammation, bruxism, inadequate oral hygiene, intra-oral infection, insufficient bone quality, peri-implant infection, endodontic or periodontal problems, untreated periodontitis, smoking and allergy to chemical elements comprising its raw material, such as V, Ti, Al, O, N, Fe, H and C.

## 8. Cleaning of the Components

Rielens makes available all of its Sanitized Abutments and Mini Abutments and can proceed to the autoclave sterilization process and used on the patient. If the Abutments are worked or modified, the cleaning procedure should follow the Practical Guide, on the website: [www.riellens.com.br](http://www.riellens.com.br).

Final restoration should be cleaned and disinfected as directed by the restorative material manufacturer prior to intraoral use.

In the case of returns, any used components that must be returned to the after sales service must be shipped sterile after pre-disinfection, cleaning and decontamination in the office.

## 9. Sterilization

The Conical Abutment, Mini Conical Abutment, Interface, Tibase Abutment, Temporary Abutment, Immediate Temporary Abutment, Abutment, Solid Abutment and Screws are non-sterile, so they should be sterilized according to professional protocol and manufacturer's instructions prior to use used on the patient. Rielens recommends placing the components one by one in a self-adhesive envelope for autoclaving steam sterilization, wet heat at 134 °C +/- 1 °C for 18 minutes, and reesterilization is prohibited, may cause maladaptation and peri-implantitis

**ATTENTION**  
Dental Product Intended for the use of qualified professionals. Supplied **NON-STERILE** and must be sterilized before use.

## 10. Precautions and Warnings

Failure to follow the procedures described in this instruction for use may lead to one or all of the following complications: implant or component damage; loosening of components; inadequate final restoration or maladaptation of mobile and fixed prostheses; implant failure. Product of SINGLE USE, reuse of components will result in patient health damage that may come into contact with a contaminated product. Use of non-compatible components or instruments will void all warranties. Use and proper handling of this product are the sole responsibility of the user. Special attention should be given to patients with chronic periodontal inflammation, bruxism, inadequate oral hygiene, intraoral infection, poor bone quality, peri-implant infection, endodontic or periodontal problems, untreated periodontitis, smoking. Rielens assumes no liability for damage resulting from the misuse of its products. The installation of the prosthesis should be done in a surgical environment and should be clean and suitable for the handling and use of health products.

## 11. Warnings

Never exceed the maximum allowable tightening torque for the abutment screw. Torison greater than the recommended torque may result in damage to the screw fit and fracture causing prosthesis failure.  
Torsion below the recommended torque may result in loosening of the component that may lead to prosthesis failure.  
Use of components that have corrosion characteristics may cause toxic effects to the patient.  
Possible complications:  
- Pain, swelling, fever, mucositis, phonetic difficulties, gingival inflammation, localized or systemic infections, aesthetic problems, restoration or dental fractures.  
- Pain, mucositis, gingival inflammation, localized or systemic infections, aesthetic problems, restoration or dental fractures.  
- Pain, swelling, fever, mucositis, phonetic difficulties, gingival inflammation, localized or systemic infections, aesthetic problems, restoration or dental fractures.

## 12. Attention

- Use the calibrated torquemeter until the recommended torque is obtained to tighten the screw and use the fixing wrench, according to table 2, to stabilize. Follow the Torquemeter manufacturer's instructions.

Ensure that the connection screws between the implant and the abutment are correctly aligned on all secondary components to prevent wear on the threads.

- Excess cement left inadvertently is a major problem and can result in soft tissue damage, bone loss and / or chronic inflammation.

- When used intra-orally, the products must be protected against any risk of aspiration. Protect components against aspiration (eg with a throat plug or floss).

## 13. Post-implementation Care

Advise the patient to see the dental surgeon periodically.

## 14. Condition of Storage

Always store at room temperature between 15°C and 35°C in the original package. Do not expose to direct sunlight.

## 15. Disposal of the Product

For disposal of Abutments and Mini Abutments used or considered inadequate, it is necessary to follow the legal sanitary procedures in force for the disposal of solid wastes in the health area and must be uncharacterized, they may be cut, bent or filed.

## 16. MRI - Magnetic Resonance Safety

The non-clinical test demonstrated that the product is MR conditional. A patient with this device may be safely checked on an MR system that meets the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla only.

- Magnetic wires of maximum spatial gradient of 4000 Gauss/cm or minimum of (40T/m).

- Maximum reported MR system, full body average specific absorption rate (SAR) of 2W/kg (Normal Operating Mode) or 4W/kg (First Level Controlled Mode).

Under the scan conditions defined above, the product is expected to produce a maximum temperature rise of 4.1°C after 15 minutes of continuous operation. In non-clinical tests, the image artifact caused by the device extends approximately 30mm from the product when represented with a sequence of gradient echo pulses with a magnetic resonance system of 3 Tesla.

If there is no MR symbol on the product label, note that the product has not been evaluated for safety and compatibility in the MR environment. The product has not been tested for heating or migration in the MR environment.

ATTENTION: The information above was developed through bibliographical equivalence of products on the market. The system has not been tested by Rielens for safety and compatibility with MRI scans. The system has not been tested for warming or migration on MRI scans. Professional tell your patient if he or she needs to perform an MRI on the above specifications.

## 17. Expiration Date and Batch

See the manufacture date, expiration and batch on the package.




## 18. INSTRUCTION FOR USE WARNING

This INSTRUCTION FOR USE is available in non-printed (electronic) format, through the electronic address of the manufacturer at [www.riellens.com.br](http://www.riellens.com.br). The INSTRUCTIONS FOR USE are indexed on the website through the respective product TRADE NAME, informed on the label of the product purchased.

If it is in the professional's interest, the INSTRUCTIONS FOR USE may be provided in printed form, at no additional cost. We ask the professionals to notify and inform Rielens of any problems that may occur with the product during its use. The requests must be made to the manufacture's SAC (Customer Service – Telephone 55 11 – 4158-5018).

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