ATENÇÃO

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WARNING

The use of Riellen's Healing and Protection Covers is specifically for use by qualified professionals with specific training in prosthetics and implantology, including diagnosis, planning and surgical protocol. They are indicated for surgical use. The Healing and Protection Covers provided by Riellen's are offered in **Non-sterile** condition.
 Never use Riellen's Healing and Protection Caps without performing the sterilization process previously.

1 – Intentions of Use

The Healing and Protective Covers have the function of preserving the platform and the inside of the implant, preserving, and increasing, if necessary, the inserted gum and adjust the thickness of soft tissues, healing the peri-implant gingival tissue.

2- Indication for Use

The indication for use of the Healing and Protection Covers is guided according to the following table:

HEALING COMPONENTS/PROTECTIVE COVER											
COMPONENTS	INDICATION FOR USE	FEATURES	CONNECTI ON	DIAMETER	HEIGHT mm						
					2	3	4	5	6	7	
	The Healer and the protective cover are used in the second phase of implant surgery, with the purpose of preserving the implant platform, creating a cylindrical space in the gingiva, carrying out adequate healing of the gingival tissue.	It is a single body component applied directly on the implant	HE	NP, RP, and WP	x	x	x	x	x	x	
			HI	NP, RP, WP, Ø3.75, 4.1, 4.3, 5.0, and 6.0	x	x	x	x	x	x	
			TI	NP, RP, WP, and 6.0	x	x	x	x	x	x	
			CM	Ø2.9, 3.3, 3.8, 4.0, 4.5, RN, WN, NC, RC, NP, WP	x	x	x	x	x	x	
		HE	NP, RP, and WP	It is a single body component applied directly on the intermediate part	NP, RP, WP, Ø3.75, and 6.0	NP, RP, WP and 6.0	2.9, RN, WN, NC, RC, NP, and WP	x	x	x	x
		HI	NP, RP, WP, Ø3.75, and 6.0								
		TI	NP, RP, WP and 6.0								
		CM	2.9, RN, WN, NC, RC, NP, and WP								

Key: HE: External Hexagon; HI: Internal Hexagon; CM: Morse Cone; TI: Inner Triangle

3 – Technical specifications and characteristics



Riellen's Healing and Protection Covers are manufactured in Titanium according to ASTM F136 and ISO 5832-3 (Standard Specification For Wrought Titanium-6Aluminium-4Vanadium). They are cylindrical in shape and have a hexagonal socket (1.2mm) or Unigrip on the upper end for the use of instruments and are provided non-sterile.

4 – Way of Presentation

Packaging

BE CAREFUL WHEN OPENING THE PACKAGE, PREVENTING FALLS, DEFORMATIONS, AND THE ACCUMULATION OF DIRT THAT MIGHT INTERFERE IN THE FINAL APPLICATION OF THE PRODUCT.

The Healing and Covers are individually packaged in PET blister packs with surgical grade paper ensuring the integrity and cleanliness until the professional uses them. Information such as description and fitting are available on the product label.

5 – Handling

To open the package, hold the blister firmly and pull the surgical paper, as indicated in the image above. The Healer and protective Cover must be used for a short term, and it is recommended to remove them within 30 days after implantation. For greater safety, protect surgical drape, avoiding that if the product is released into the mouth, it will not be swallowed or aspirated by the patient.

– Position the Healer/Protective Cover correctly over the implant connection and apply the recommended manual torque of 15Ncm, making sure that there are no micro gaps that allow the transit of microorganisms.

6 – Contraindications

The use of Riellen's Healing and protective covers are contraindicated in cases of patients who are allergic to any chemical element that make the raw material such as: V, Ti, Al, O, N, Fe, H and C.

7- Adverse Effects/Technical Complaints

Possible complications: Gingival inflammation, localized infections, maladaptation, mechanical problems

8-Component cleaning

Riellen's makes available all its Healing and Protective Covers Sanitized.

9 – Sterilization

These products are provided non-sterile; thus, they shall be sterilized in an autoclave.

- Use chemical and biological indicators for each sterilization cycle, according to practice protocol.
- Allow Riellen's Healing /Protective Covers to dry and cool in the autoclave before handling to avoid contamination.
- Riellen's recommends sterilization at the parameters of 134 °C for 15 minutes.
- Riellen's recommends placing the components in their own packaging for sterilization in an autoclave.

10 – Precautions and Warnings

– The use of components or instruments that are not compatible will render all warranties void. The correct use and handling of this product is the sole responsibility of the user. Riellen's takes no responsibility for damages resulting from incorrect use of its products.

-The environment for installing the prosthesis must be clean and suitable for the handling and use of health products.

– Healing and Protective covers are instruments of the **"FORBIDDEN TO REPROCESS"** type; the reuse of the components will result in damage to the patient's health who may come in contact with a contaminated product.

– Torque superior to the recommendation may result in damage to the socket, fracture or locking of the Heaters and Protection Covers, causing implant failure. Torque lower than the recommendation may result in component loosening, which may lead to poor healing of the gingival tissue.

WARNING

Dental product intended for use by qualified professionals. Provided NON-STERILE, never reuse Riellen's Healing /Protective Cover, risk of infection with infectious agents. Caution when Handling the component, risk of accidentally aspirating or swallowing.

11 – Alerts/Cautions

In case of technical complaints or adverse events with Riellen's products, please notify us through SAC (Customer Service).
– Ensure that the Healing and protective covers are correctly aligned to prevent thread wear.
– In the event of damage to the package (blister) or after the expiration date, Riellen's Healing and Protective covers cannot be used.
– Any used components that must be returned to the after-sales service should be shipped sterile after pre-disinfection, cleaning, and decontamination in the practice.

12 – Storage Conditions

Always store at room temperature between 15 °C and 35 °C in the original package. Do not expose to direct sunlight.

13 – Product Disposal

For the disposal of the Healing and Protective covers used or considered inadequate, it is necessary to follow the current legal procedures for the disposal of products, and they must be defused by being cut, bent, or filed away.


14 – Expiration Date and Batch

See the manufacture date, expiration, and batch on the package.

15 – ALERT INSTRUCTION FOR USE

This INSTRUCTION FOR USE is available in non-printed (electronic) format, through the electronic address of the manufacturer at <http://www.riellens.com.br/>. The INSTRUCTIONS FOR USE are indexed on the website through the respective product TRADE NAME, informed on the label of the product purchased. If it is in the professional's interest, the INSTRUCTIONS FOR USE may be provided in printed form, at no additional cost. The request must be made to the SAC (Customer Service) of the manufacturer.

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Technician in Charge/Technical Responsible:
Vanessa Guimarães Souza de Castro
CRFSP No. / Professional ID: 63664
ANVISA Registration No.: 80523780005

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