ATENÇÃO

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WARNING



Riellen's Transfers are specifically for use by qualified professionals with specific training in prosthetics and implantology, including diagnosis, planning and surgical protocol. They are indicated for laboratory and surgical use. The Transfers provided by Riellen's are offered in **non-sterile condition**. Never use Riellen's Transfers without performing the sterilization process.

1 – Intentions of Use

The function of the Riellen's Transfer is to transfer the position of the implant and intermediates (Conical Abutment / Mini Abutment / Universal Abutment) so that the prosthesis can be manufactured. It is indicated for open tray (drag) or closed tray (replacement) techniques.

2 – Indication for Use

The indication of the use of the Transfer is guided according to the following table:

RIELLEN'S TRANSFER COMPONENTS					
COMPONENT	INDICATION FOR USE	Models	FEATURES	CONNECTION	DIAMETER
 Stainless Steel Transfer	Used to transfer the positioning of the dental implant installed in the oral cavity to the plaster model, by means of an impression.	Open Tray	Used for impressions of single and multiple prosthetics with open tray	HE	NP, RP, and WP
				HI	Ø3.3 and 3.75
		Closed Tray	Used for impressions of single and multiple prosthetics with closed tray	TI	Ø3.3
				CM	Ø3.0, 3.8, 4.5, Ø4.5/5.0 RN, and WN
 Titanium Transfer		Open Tray	Used for impressions of single and multiple prosthetics with open tray	HI	NP, RP, and WP
				TI	Ø3.5/4.0, 4.5/5.0, NP, RP, and WP
		Closed tray	Used for impressions of single and multiple prosthetics with closed tray	CM	Ø2.9, 3.8, NN, NP, RN, RP, WN, NC and RC

Key: HE: External Hexagon; HI: Internal Hexagon; CM: Morse Cone; TI: Inner Triangle

3 – Technical specifications and characteristics



Riellen's transfers are made of titanium according to standards ASTM F136 and ISO 5832-3 (Standard Specification for Wrought Titanium-6Aluminum-4Vanadium) or Stainless-Steel ASTM F899-09 303 (Standard Specification for Wrought Stainless Steels for Surgical Instruments) ISO 5832-1 (Standard Specification for Implants for Surgery – Metallic materials Wrought Stainless Steels), intended for single use, its reuse is forbidden, provided non-sterile and have a cylindrical in shape. They are presented in two parts, transfer + Screw and are available as open tray and closed tray in relation to the impression technique.

4 – Way of Presentation

Packaging

BE CAREFUL WHEN OPENING THE PACKAGE, PREVENTING FALLS, DEFORMATIONS, AND THE ACCUMULATION OF DIRT THAT MIGHT INTERFERE IN THE FINAL APPLICATION OF THE PRODUCT.

Riellen's Transfers are individually packaged in blister packaging ensuring the integrity and cleanliness on the product label until the professional uses them. Information such as description, size, manufacture date and batch are available on the product label.

5- Handling

To open the package, hold the blister firmly and pull the surgical paper, as indicated in the image above. Place the selected Transfer according to the planned impression technique and install directly on the implant or intermediate. To tighten the screw, use manual torque with digital wrench. Use the technique of impression of the patient's mouth and pour the plaster to obtain the replica of the implant in the oral cavity.

6- Contraindications

This product has no contraindications, as they are used correctly for the indicated purposes.

7-Component cleaning

Riellen's makes available all its Transfers Sanitized.

8- Sterilization

These products are provided non-sterile; thus, they shall be sterilized in an autoclave.

- Use chemical and biological indicators for each sterilization cycle, according to practice protocol.
- Allow Riellen's Transfers to dry and cool in the autoclave before handling to avoid contamination.
- Riellen's recommends sterilization at the parameters of 134 °C for 15 minutes.
- Riellen's recommends placing the components in their own packaging for sterilization in an autoclave.

9 – Precautions and Warnings

The Use of components or instruments that are not compatible will render all warranties void. – The correct use and handling of this product is the sole responsibility of the user. Riellen's takes no responsibility for damages resulting from incorrect use of its products. The Environment for installing the prosthesis must be clean and suitable for the handling and use of health products.

WARNING

Dental product intended for use by qualified professionals. Provided NON-STERILE, never reuse Riellen's Transfers, risk of infection with infectious agents.

10 – Caution

To tighten the screw, use a digital wrench applying the necessary torque for a fair tightening of the screw without damaging it. Ensure that the Transfers screws are correctly aligned between the implant and/or intermediate to avoid thread wear.

11 – Adverse Effects

This product has no adverse effects, if used correctly for the indicated purposes.

12 – Storage Conditions

Always store in the original package, do not expose to direct sunlight and store at a temperature of 15 °C – 35 °C.

13 – Disposal of the product

For the disposal of Transfers used or considered inadequate, it is necessary to follow the current legal procedures for the disposal of products, and they must be defused by being cut, bent, or filed away.

14 – Expiration Date and Batch

See the manufacture date, expiration, and batch on the package. This product shall be disposed if it loses its functionality.

15 – ALERT INSTRUCTION FOR USE

This INSTRUCTION FOR USE is available in non-printed (electronic) format, through the electronic address of the manufacturer at www.riellens.com.br. The INSTRUCTIONS FOR USE are indexed on the website through the respective product TRADE NAME, informed on the label of the product purchased. If it is in the professional's interest, the INSTRUCTIONS FOR USE may be provided in printed form, at no additional cost. The request must be made to the SAC (Customer Service) of the manufacturer.

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Technician In Charge/ Technical Responsible:
Vanessa Guimarães Souza de Castro
CRFSP No. / Professional ID: 63664
Anvisa Registration: Transfer No. 80523780002 and No. 805237800

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